

Thymidine phosphorylase

Catalog No: #AP89082



Package Size: #AP89082-1 50ug #AP89082-2 100ug #AP89082-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Thymidine phosphorylase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | TP, Gliostatin, Platelet-derived endothelial cell growth factor, PD-ECGF, TdRPase, ECGF1 |
| Accession No. | Uniprot ID:P19971 |
| Uniprot | P19971 |
| GeneID | 1890; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TYMP

Note: This product is for in vitro research use only