

## Transcriptional coactivator YAP1

Catalog No: #AP89084



Package Size: #AP89084-1 50ug #AP89084-2 100ug #AP89084-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Transcriptional coactivator YAP1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Yes-associated protein 1,65 kDa Yes-associated protein,YAP65,Protein yorkie homolog,YAP65
Accession No.	Uniprot ID:P46936
Uniprot	P46936
GeneID	396171;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:YAP1

Note: This product is for in vitro research use only