

## Zona pellucida sperm-binding protein 3

Catalog No: #AP89088



Package Size: #AP89088-1 50ug #AP89088-2 100ug #AP89088-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Zona pellucida sperm-binding protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Sperm receptor,Zona pellucida glycoprotein 3,Zp-3,Zona pellucida protein C,Zp-3,Zpc
Accession No.	Uniprot ID:P10761
Uniprot	P10761
GeneID	22788;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Zp3

Note: This product is for in vitro research use only