

WAP four-disulfide core domain protein 2

Catalog No: #AP89090



Package Size: #AP89090-1 50ug #AP89090-2 100ug #AP89090-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	WAP four-disulfide core domain protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Epididymal secretory protein E4,CE4,Major epididymis-specific protein E4
Accession No.	Uniprot ID:Q28894
Uniprot	Q28894
GeneID	403919;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:WFDC2

Note: This product is for in vitro research use only