

Troponin I, fast skeletal muscle

Catalog No: #AP89097



Package Size: #AP89097-1 50ug #AP89097-2 100ug #AP89097-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Troponin I, fast skeletal muscle
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Troponin I, fast-twitch isoform
Accession No.	Uniprot ID:P68246
Uniprot	P68246
GeneID	396386;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TNNI2

Note: This product is for in vitro research use only