

# Tumor necrosis factor receptor superfamily member 16

Catalog No: #AP89098

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP89098-1 50ug #AP89098-2 100ug #AP89098-3 1mg

## Description

Product Name	Tumor necrosis factor receptor superfamily member 16
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Low affinity neurotrophin receptor p75NTR,Low-affinity nerve growth factor receptor,NGF receptor,CD271,Tnfrsf16
Accession No.	Uniprot ID:Q9Z0W1
Uniprot	Q9Z0W1
GeneID	18053;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ngfr

Note: This product is for in vitro research use only