

## Tumor necrosis factor receptor superfamily member 4

Catalog No: #AP89099



Package Size: #AP89099-1 50ug #AP89099-2 100ug #AP89099-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tumor necrosis factor receptor superfamily member 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ACT35 antigen,OX40L receptor,TAX transcriptionally-activated glycoprotein 1 receptor,CD134,TXGP1L
Accession No.	Uniprot ID:P43489
Uniprot	P43489
GeneID	7293;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:TNFRSF4

Note: This product is for in vitro research use only