

## Vitamin D3 receptor

Catalog No: #AP89102



Package Size: #AP89102-1 50ug #AP89102-2 100ug #AP89102-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Vitamin D3 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	VDR,1,25-dihydroxyvitamin D3 receptor,Nuclear receptor subfamily 1 group I member 1,NR111
Accession No.	Uniprot ID:A3RGC1
Uniprot	A3RGC1
GeneID	396628;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:VDR

Note: This product is for in vitro research use only