

DNA topoisomerase 2-beta

Catalog No: #AP89103



Package Size: #AP89103-1 50ug #AP89103-2 100ug #AP89103-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	DNA topoisomerase 2-beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	DNA topoisomerase II, beta isozyme
Accession No.	Uniprot ID:Q64511
Uniprot	Q64511
GeneID	21974;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Top2b

Note: This product is for in vitro research use only