

Prosalusin

Catalog No: #AP89104



Package Size: #AP89104-1 50ug #AP89104-2 100ug #AP89104-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Prosalusin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Torsin family 2 member A,Torsin-2A
Accession No.	Uniprot ID:P0C7W2
Uniprot	P0C7W2
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tor2a

Note: This product is for in vitro research use only