

Uveal autoantigen with coiled-coil domains and ankyrin repeats

Catalog No: #AP89108

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89108-1 50ug #AP89108-2 100ug #AP89108-3 1mg

Description

Product Name	Uveal autoantigen with coiled-coil domains and ankyrin repeats
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Nuclear membrane-binding protein, Nucling, Kiaa1561
Accession No.	Uniprot ID: Q8CGB3
Uniprot	Q8CGB3
GeneID	72565;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Uaca

Note: This product is for in vitro research use only