

UDP-glucuronosyltransferase 1-1

Catalog No: #AP89111



Package Size: #AP89111-1 50ug #AP89111-2 100ug #AP89111-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	UDP-glucuronosyltransferase 1-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	UDPGT 1-1,B1,UDP-glucuronosyltransferase 1A1,Ugt1
Accession No.	Uniprot ID:Q64550
Uniprot	Q64550
GeneID	24861;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ugt1a1

Note: This product is for in vitro research use only