

Sortilin

Catalog No: #AP89113



Package Size: #AP89113-1 50ug #AP89113-2 100ug #AP89113-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sortilin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Glycoprotein 110,Gp110,Neurotensin receptor 3,NTR3
Accession No.	Uniprot ID:O54861
Uniprot	O54861
GeneID	83576;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sort1

Note: This product is for in vitro research use only