

## Urocortin-3

Catalog No: #AP89114



Package Size: #AP89114-1 50ug #AP89114-2 100ug #AP89114-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Urocortin-3  |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Ms   |
| Other Names        | Urocortin III,Ucn III                                  |
| Accession No.      | Uniprot ID:Q924A4                                      |
| Uniprot            | Q924A4   |
| GeneID             | 83428;   |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ucn3

Note: This product is for in vitro research use only