

Mitochondrial uncoupling protein 3

Catalog No: #AP89115



Package Size: #AP89115-1 50ug #AP89115-2 100ug #AP89115-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitochondrial uncoupling protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	UCP 3,Solute carrier family 25 member 9,SLC25A9
Accession No.	Uniprot ID:O77792
Uniprot	O77792
GeneID	281563;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:UCP3

Note: This product is for in vitro research use only