

Vascular endothelial growth factor B

Catalog No: #AP89116



Package Size: #AP89116-1 50ug #AP89116-2 100ug #AP89116-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vascular endothelial growth factor B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	VEGF-B,VEGF-related factor,VRF,VRF
Accession No.	Uniprot ID:P49765
Uniprot	P49765
GeneID	7423;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:VEGFB

Note: This product is for in vitro research use only