

TYRO protein tyrosine kinase-binding protein

Catalog No: #AP89120



Package Size: #AP89120-1 50ug #AP89120-2 100ug #AP89120-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	TYRO protein tyrosine kinase-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	DNAX-activation protein 12,DAP12
Accession No.	Uniprot ID:Q95J79
Uniprot	Q95J79
GeneID	282390;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TYROBP

Note: This product is for in vitro research use only