

Vigilin

Catalog No: #AP89122



Package Size: #AP89122-1 50ug #AP89122-2 100ug #AP89122-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Vigilin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	High density lipoprotein-binding protein,HDL-binding protein,Hbp
Accession No.	Uniprot ID:Q9Z1A6
Uniprot	Q9Z1A6
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hdlbp

Note: This product is for in vitro research use only