

WAP, Kazal, immunoglobulin, Kunitz and NTR domain-containing protein 1

Catalog No: #AP89126

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89126-1 50ug #AP89126-2 100ug #AP89126-3 1mg

Description

Product Name	WAP, Kazal, immunoglobulin, Kunitz and NTR domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Growth and differentiation factor-associated serum protein 2,GASP-2,WAP, follistatin, immunoglobulin, Kunitz and NTR domain-containing protein,C16orf12,GASP2,WFIKKN,PSEC0040
Accession No.	Uniprot ID:Q96NZ8
Uniprot	Q96NZ8
GeneID	117166;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:WFIKKN1

Note: This product is for in vitro research use only