

Tumor necrosis factor receptor superfamily member 8

Catalog No: #AP89130



Package Size: #AP89130-1 50ug #AP89130-2 100ug #AP89130-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CD30L receptor, Lymphocyte activation antigen CD30, CD30, Cd30
Accession No.	Uniprot ID:P97525
Uniprot	P97525
GeneID	25069;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tnfrsf8

Note: This product is for in vitro research use only