

Serotransferrin

Catalog No: #AP89135



Package Size: #AP89135-1 50ug #AP89135-2 100ug #AP89135-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Serotransferrin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Transferrin,Beta-1 metal-binding globulin,Liver regeneration-related protein LRRG03,Siderophilin
Accession No.	Uniprot ID:P12346
Uniprot	P12346
GeneID	24825;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tf

Note: This product is for in vitro research use only