

WNT1-inducible-signaling pathway protein 1

Catalog No: #AP89136



Package Size: #AP89136-1 50ug #AP89136-2 100ug #AP89136-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	WNT1-inducible-signaling pathway protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	WISP-1,CCN family member 4,ELM-1,Ccn4,Elm1
Accession No.	Uniprot ID:Q99PP0
Uniprot	Q99PP0
GeneID	65154;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Wisp1

Note: This product is for in vitro research use only