

Ubiquitin-associated protein 2

Catalog No: #AP89140



Package Size: #AP89140-1 50ug #AP89140-2 100ug #AP89140-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ubiquitin-associated protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	UBAP-2,KIAA1491
Accession No.	Uniprot ID:Q5T6F2
Uniprot	Q5T6F2
GeneID	55833;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:UBAP2

Note: This product is for in vitro research use only