

Urea transporter 1

Catalog No: #AP89146



Package Size: #AP89146-1 50ug #AP89146-2 100ug #AP89146-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Urea transporter 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Solute carrier family 14 member 1,Urea transporter B,UT-B,Urea transporter, erythrocyte,UT11,UT3
Accession No.	Uniprot ID:P97689
Uniprot	P97689
GeneID	54301;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc14a1

Note: This product is for in vitro research use only