

V-type proton ATPase catalytic subunit A

Catalog No: #AP89147



Package Size: #AP89147-1 50ug #AP89147-2 100ug #AP89147-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	V-type proton ATPase catalytic subunit A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	V-ATPase subunit A,V-ATPase 69 kDa subunit,Vacuolar proton pump subunit alpha
Accession No.	Uniprot ID:Q90647
Uniprot	Q90647
GeneID	395821;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ATP6V1A

Note: This product is for in vitro research use only