

Troponin C, slow skeletal and cardiac muscles

Catalog No: #AP89148



Package Size: #AP89148-1 50ug #AP89148-2 100ug #AP89148-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Troponin C, slow skeletal and cardiac muscles
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	TN-C, TNNC
Accession No.	Uniprot ID:P63315
Uniprot	P63315
GeneID	509486;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TNNC1

Note: This product is for in vitro research use only