

Troponin T, slow skeletal muscle

Catalog No: #AP89149



Package Size: #AP89149-1 50ug #AP89149-2 100ug #AP89149-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Troponin T, slow skeletal muscle
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	TnTs, Slow skeletal muscle troponin T, sTnT
Accession No.	Uniprot ID:Q75ZZ6
Uniprot	Q75ZZ6
GeneID	396579;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TNNT1

Note: This product is for in vitro research use only