

OR52B6 Antibody

Catalog No: #44587

Package Size: #44587-1 50ul #44587-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	OR52B6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Human;Rat;Mouse;
Specificity	OR52B6 Antibody detects endogenous levels of total OR52B6
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human OR52B6
Target Name	OR52B6
Other Names	OR52B6, OR11-47, Olfactory receptor 52B6, Olfactory receptor OR11-47
Accession No.	Swiss-Prot#: Q8NGF0NCBI Gene ID:340980
Uniprot	Q8NGF0
GeneID	340980;
Calculated MW	37 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only