

Vascular endothelial growth factor C

Catalog No: #AP89159



Package Size: #AP89159-1 50ug #AP89159-2 100ug #AP89159-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vascular endothelial growth factor C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	VEGF-C,Flt4 ligand,Flt4-L,Vascular endothelial growth factor-related protein,VRP
Accession No.	Uniprot ID:O35757
Uniprot	O35757
GeneID	114111;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Vegfc

Note: This product is for in vitro research use only