

## Vacuole membrane protein 1

Catalog No: #AP89161



Package Size: #AP89161-1 50ug #AP89161-2 100ug #AP89161-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Vacuole membrane protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	NF-E2-inducible protein 2,Protein ni-2,Transmembrane protein 49,Tmem49
Accession No.	Uniprot ID:Q99KU0
Uniprot	Q99KU0
GeneID	75909;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Vmp1

Note: This product is for in vitro research use only