

## AN1-type zinc finger protein 6

Catalog No: #AP89162



Package Size: #AP89162-1 50ug #AP89162-2 100ug #AP89162-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	AN1-type zinc finger protein 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Associated with PRK1 protein,Zinc finger A20 domain-containing protein 3,AWP1,ZA20D3,HT032
Accession No.	Uniprot ID:Q6FIF0
Uniprot	Q6FIF0
GeneID	54469;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ZFAND6

Note: This product is for in vitro research use only