

## Tumor necrosis factor ligand superfamily member 13B

Catalog No: #AP89168



Package Size: #AP89168-1 50ug #AP89168-2 100ug #AP89168-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tumor necrosis factor ligand superfamily member 13B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	B lymphocyte stimulator,BLYS,B-cell-activating factor,BAFF,Dendritic cell-derived TNF-like molecule,TNF- and APOL-related leukocyte expressed ligand 1,TALL-1,CD257,BAFF,BLYS,TALL1,TNFSF20,ZTNF4,UNQ401/PRO738
Accession No.	Uniprot ID:Q9Y275
Uniprot	Q9Y275
GeneID	10673;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:TNFSF13B

Note: This product is for in vitro research use only