

## Tumor necrosis factor ligand superfamily member 18

Catalog No: #AP89169



Package Size: #AP89169-1 50ug #AP89169-2 100ug #AP89169-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tumor necrosis factor ligand superfamily member 18
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	GITR ligand,GITRL,Glucocorticoid-induced TNF-related ligand,Gitrl
Accession No.	Uniprot ID:Q7TS55
Uniprot	Q7TS55
GeneID	240873;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Tnfsf18

Note: This product is for in vitro research use only