

Tumor necrosis factor

Catalog No: #AP89170



Package Size: #AP89170-1 50ug #AP89170-2 100ug #AP89170-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cachectin, TNF-alpha, Tumor necrosis factor ligand superfamily member 2, TNF-a, Tnfa, Tnfsf2
Accession No.	Uniprot ID:P06804
Uniprot	P06804
GeneID	21926;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tnf

Note: This product is for in vitro research use only