

Tumor necrosis factor ligand superfamily member 6

Catalog No: #AP89173



Package Size: #AP89173-1 50ug #AP89173-2 100ug #AP89173-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor ligand superfamily member 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	CD95 ligand,CD95-L,Fas antigen ligand,Fas ligand,CD178,CD95L,FASL,TNFSF6
Accession No.	Uniprot ID:Q9BEA8
Uniprot	Q9BEA8
GeneID	396726;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FASLG

Note: This product is for in vitro research use only