

Thioredoxin reductase 1, cytoplasmic

Catalog No: #AP89176



Package Size: #AP89176-1 50ug #AP89176-2 100ug #AP89176-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Thioredoxin reductase 1, cytoplasmic
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	TR,NADPH-dependent thioredoxin reductase,Thioredoxin reductase TR1,Trxr1
Accession No.	Uniprot ID:O89049
Uniprot	O89049
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Txnrd1

Note: This product is for in vitro research use only