

Tryptase beta-2

Catalog No: #AP89178



Package Size: #AP89178-1 50ug #AP89178-2 100ug #AP89178-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tryptase beta-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Tryptase-2, Mast cell protease 6, mMCP-6, Mcpt6
Accession No.	Uniprot ID:P21845
Uniprot	P21845
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tpsb2

Note: This product is for in vitro research use only