

Transcription termination factor 2

Catalog No: #AP89179



Package Size: #AP89179-1 50ug #AP89179-2 100ug #AP89179-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcription termination factor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Lodestar homolog, RNA polymerase II termination factor, Transcription release factor 2, F2
Accession No.	Uniprot ID:Q9UNY4
Uniprot	Q9UNY4
GeneID	8458;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TTF2

Note: This product is for in vitro research use only