

WNT1-inducible-signaling pathway protein 2

Catalog No: #AP89180



Package Size: #AP89180-1 50ug #AP89180-2 100ug #AP89180-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	WNT1-inducible-signaling pathway protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	WISP-2,CCN family member 5,CCN family protein COP-1,Connective tissue growth factor-like protein,CTGF-L,Ccn5,Cop1,Ctgfl
Accession No.	Uniprot ID:Q9JHC6
Uniprot	Q9JHC6
GeneID	29576;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Wisp2

Note: This product is for in vitro research use only