

Protein Wnt-5b

Catalog No: #AP89181



Package Size: #AP89181-1 50ug #AP89181-2 100ug #AP89181-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein Wnt-5b
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Wnt-5b
Accession No.	Uniprot ID:P22726
Uniprot	P22726
GeneID	22419;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Wnt5b

Note: This product is for in vitro research use only