

Somatostatin receptor type 2

Catalog No: #AP89185



Package Size: #AP89185-1 50ug #AP89185-2 100ug #AP89185-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Somatostatin receptor type 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	SS-2-R,SS2-R,SS2R
Accession No.	Uniprot ID:Q49LX6
Uniprot	Q49LX6
GeneID	403472;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SSTR2

Note: This product is for in vitro research use only