

Zinc-alpha-2-glycoprotein

Catalog No: #AP89187



Package Size: #AP89187-1 50ug #AP89187-2 100ug #AP89187-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Zinc-alpha-2-glycoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Zn-alpha-2-GP
Accession No.	Uniprot ID:Q3ZCH5
Uniprot	Q3ZCH5
GeneID	508800;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: AZGP1

Note: This product is for in vitro research use only