

Tumor necrosis factor ligand superfamily member 12

Catalog No: #AP89188



Package Size: #AP89188-1 50ug #AP89188-2 100ug #AP89188-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor ligand superfamily member 12
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	TNF-related weak inducer of apoptosis,TWEAK
Accession No.	Uniprot ID:O54907
Uniprot	O54907
GeneID	21944;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tnfsf12

Note: This product is for in vitro research use only