

Lymphotactin

Catalog No: #AP89192



Package Size: #AP89192-1 50ug #AP89192-2 100ug #AP89192-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lymphotactin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ATAC,C motif chemokine 1,Cytokine SCM-1,Lymphotaxin,SCM-1-alpha,Small-inducible cytokine C1,XC chemokine ligand 1,LTN,SCYC1
Accession No.	Uniprot ID:P47992
Uniprot	P47992
GeneID	6375;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:XCL1

Note: This product is for in vitro research use only