

Urokinase plasminogen activator surface receptor

Catalog No: #AP89196



Package Size: #AP89196-1 50ug #AP89196-2 100ug #AP89196-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Urokinase plasminogen activator surface receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	U-PAR, Monocyte activation antigen Mo3, CD87, MO3, UPAR
Accession No.	Uniprot ID: Q03405
Uniprot	Q03405
GeneID	5329;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PLAUR

Note: This product is for in vitro research use only