

Zinc finger E-box-binding homeobox 1

Catalog No: #AP89197



Package Size: #AP89197-1 50ug #AP89197-2 100ug #AP89197-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Zinc finger E-box-binding homeobox 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Transcription factor 8,TCF-8,Zinc finger homeodomain enhancer-binding protein,Zfhpf,Tcf8
Accession No.	Uniprot ID:Q62947
Uniprot	Q62947
GeneID	25705;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Zeb1

Note: This product is for in vitro research use only