

Pantetheinase

Catalog No: #AP89198



Package Size: #AP89198-1 50ug #AP89198-2 100ug #AP89198-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pantetheinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Pantetheine hydrolase,Vascular non-inflammatory molecule 1, Vanin-1
Accession No.	Uniprot ID:Q58CQ9
Uniprot	Q58CQ9
GeneID	526704;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:VNN1

Note: This product is for in vitro research use only