

## WAP four-disulfide core domain protein 2

Catalog No: #AP89200

Package Size: #AP89200-1 50ug #AP89200-2 100ug #AP89200-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	WAP four-disulfide core domain protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Epididymal secretory protein E4
Accession No.	Uniprot ID:Q8MI69
Uniprot	Q8MI69
GeneID	396987;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: WFDC2

Note: This product is for in vitro research use only