

Protein Wnt-2b

Catalog No: #AP89202



Package Size: #AP89202-1 50ug #AP89202-2 100ug #AP89202-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein Wnt-2b
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Protein Wnt-13,WNT13
Accession No.	Uniprot ID:Q93097
Uniprot	Q93097
GeneID	7482;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:WNT2B

Note: This product is for in vitro research use only