

E3 ubiquitin-protein ligase TRIM63

Catalog No: #AP89204



Package Size: #AP89204-1 50ug #AP89204-2 100ug #AP89204-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | E3 ubiquitin-protein ligase TRIM63 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | Muscle-specific RING finger protein 1,MuRF-1,RING-type E3 ubiquitin transferase TRIM63,Tripartite motif-containing protein 63,Murf1 |
| Accession No. | Uniprot ID:Q91Z63 |
| Uniprot | Q91Z63 |
| GeneID | 140939; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Trim63

Note: This product is for in vitro research use only